



# Film Reticles




Pyser-SGI film scales and grids are being used in so many different ways, far beyond their original intention :


1. As a coverslip, so giving direct measurements.
2. As a counting grid.
3. As an eyepiece reticle – just cut to the diameter required to fit in your eyepiece – guide circles are included (FNE1 only).
4. As a stage micrometer – to calibrate your microscope.
5. As a specimen support.
6. As a direct reading scale with magnifiers.




Pattern	Description	Order Code	Price
FNE1	Plastic reticle strips, 5 images per strip, 2 strips per pack - total 10 images, 150mm x 30mm x 0.18mm thick film, guide circles of 16mm, 19mm, 21mm and 24mm around each scale image for cutting	01AFS201	€47.00 per pack
CMC34A	10 x 10 grid of 0.1mm squares at 5 positions identified as A - E, indexed 1 - 10 along top and A - J down the side, 25pcs per box, each 22mm x 22mm.	01ACMC34A	€92.00 per box
CMC35	10 x 10 grid of 1mm squares, each square individually identified 0 to 99, 25pcs per box, each 22mm x 22mm.	01ACMC35	€81.00 per box
CMC71	20 x 20 grid of 0.5mm squares, indexed 1 - 20 along top, A - T down the side and on central cross, 25pcs per box, each 22mm x 22mm.	01ACMC71	€81.00 per box
CMC01	10mm scale in 0.1mm divisions, 25pcs per box, each 22mm x 22mm.	01ACMC01	€81.00 per box




CMC34A



CMC35



CMC71



CMC01